

TX10 Medserver



TX10 Medserver Features & Benefits

- Stable, trim-line footprint for simple maneuverability
- Keyless access with Auto-Locking
- Configurable medication storage with 2 bin sizes and 3 drawer sizes
- Open platform for a host of computing options including All-In-One PC, laptop or thin-client solutions
- Extended power system for run times of 12+ hours
- Full selection of convenience accessories and options
- Ergonomic design and PowerLift height adjustment
- Expansive work area with integrated forward-facing handle

TX10 Medserver Specifications

Base Dimensions:	16" x 19.5"
Work Surface Area:	18" x 13" with integrated forward-facing handles
Work Surface Height:	30" - 42"
Medserver Height (including Monitor):	49" - 61"
Construction:	Formed steel, aluminum, heavy-duty polymer work surface
Finish:	Durable powder coat
Wheels:	2 locking, 2 standard, tracking optional

Specifications, availability & components are subject to change without notice.

Contact: Christopher Sunda
Healthcare Technology Consultant
csunda@healthcarecart.com

201-406-4797 PHONE
888-573-6074 FAX

Computing Solutions



Top Mount

Available on all Mobility Platforms



Rear Mount

Available on all Mobility Platforms

Laptop PC Mounting Solutions TX10 Medserver

Artromick mobility platforms accommodate a host of laptop models from the leading healthcare IT brands including **HP, Dell, IBM, Panasonic and Toshiba.**

- Flexible mounting options
- Secure and stable mounting platforms
- Can be adjusted for different laptop sizes to permit hardware replacements and upgrades
- Available with choice of 17" or 19" Planar LCD monitor featuring vivid user interface
- Organized cable management with accommodation for bar code scanner and external wireless devices

Artromick mobile computing platforms accommodate a wide range of software and hardware solutions.

Software applications

Epic	GE Healthcare	McKesson	Eclipsys
Cerner	Keane	Misys	QuadraMed
Meditech	Picis	Siemens	InterSystems
Mediserve	NextGen	Allscripts	Wellogix

Consult your Healthcarecart.com representative for application not listed above.

Hardware systems

Dell	Panasonic	Neoware	Fujitsu
Sony	Hewlett-Packard	Gateway	Compact
IBM/Lenovo	Motion Computing	Sharp	Datalux
Toshiba	Pelham Sloane	Wyse	NEC

Contact: Christopher Sunda
Healthcare Technology Consultant
csunda@healthcarecart.com

201-406-4797 PHONE
888-573-6074 FAX

Computing Solutions



All-in-One

All-In-One Computers - TX10 Medserver

The NEW All-in-One computing solution offers high-performance processing, extended battery run times, and a host of features that promote improved infection control. Powerful and durable, Artromick All-in-One computers are ideal for supporting the latest point-of-care workflow solutions.

- Available on all Mobility Platforms
- Flexible mounting options
- Powerful Intel® Core 2 Duo processing
- Vibrant 17" screen
- Integrated wireless for reliable connectivity
- Touch screen monitor optional
- Feature sets designed specifically for the healthcare environment



Thin Client Options

Thin-Client Computing

Select from the latest thin-client solutions including Neoware, HP and WYSE. Inquire with your Artromick representative about other thin-client devices.

- Available on all Mobility Platforms
- Flexible mounting options
- Choice of monitor options
- Optional touch screens
- Supports wireless data environment
- State-of-the-art computing performance



Tablet Options

Tablet Computing

All mobility platforms accommodate tablet solutions including Dell, Toshiba & Motion Computing and other popular healthcare brands. Talk with your Sales Representative to get the best fit solution for your tablet PC.

- Flexible mounting options including rear-mount
- Supports wireless data environment
- State-of-the-art computing performance

Accessories & Options


Convenience Accessories - TX10 & TX11 Medserver




**Locking & Tracking
Casters**



Sharps Container



Storage Basket



Touchscreen Stylus




Waste Container



Mouse Holder



Sealed Keyboard



Keyboard with Cover




Slide-Out Mouse Tray



Sealed Mouse



Optical Mouse



Monitor Security Screen

Contact: Christopher Sunda
Healthcare Technology Consultant
csunda@healthcarecart.com

201-406-4797 PHONE
888-573-6074 FAX

Accessories & Options

Convenience Accessories - TX10 & TX11 Medserver



**Medserver Bin
Dividers**



Cup Holder



Monitor Pivot Mount




Scanner Holder



Glove Box



I.V. Pole



Wall Lock



Organizer Trays



Clear Top Mat

Security Options

Lock Systems & Options - TX10 & TX11 Medserver

Keyless Access with Auto-Relocking



Artromick has fitted the TX10 & TX11 Medservers with its proven and innovative automatic relocking system and the latest integrated medication cart security. To further optimize security,

add a secondary access requirement such as ID badge card, bar code or proximity reader.

Standard Features

AutoLock™

Automatically secures the cart whenever the cart remains motionless for a user-defined relock time interval of 1 to 99 minutes.

Easy-to-Edit Settings

Display screen shows messages, date and time settings, and easy-to-follow user instructions.

Multiple User Access

Users can be assigned different access rights based on job function and authorization.

Override Code

Provides authorized users with a special means of access that bypasses alarms, timers, and relocking mechanisms. Typically used while the cart is being restocked by the pharmacy.

LCD Battery Status Indicator

The extended-life battery system ensures an extra-long cycle of useful power supply with battery charge status clearly indicated in percentage on the LCD display.

Security Upgrades




Secondary Access Security

Select from a bar code, mag stripe or proximity card reader to serve as the primary or secondary access requirement for the TX10 Medserver.



Audit Tracking

Authorized users can track cart activity, including user access and relocking for over 2,000 events.



Remote Management

Wirelessly enables authorized users to program all carts and retrieve event logs, add/delete user profiles and adjust cart settings from a single PC.